MICROSCOPE

Instruction and Activity Guide

Introduce you to the hidden wonders of the microscopic world!

CAUTION: To avoid eye injury, do not reflect direct sunlight through the microscope.

Getting acquainted with the microscope

Before students begin observing specimens, have them take some time to learn about the different parts of the microscope. Their experiments will be more successful!

To begin: carefully remove the microscope from the package using both hands to grasp the back and the base of the unit. Place the microscope on a flat surface near a light source such as a lamp or a window. Next, introduce students to the microscope by using this diagram to locate the different parts:

- A. Eyepiece: has a 10X-20X zoom lens.
- B. Body Tube: extends below the eyepiece and is used to focus the lenses.
- C. Focus Knobs: are located at the arm top. Turn these two knobs and see the focus change in the eyepiece. Turn slowly to get the best focus.
- D. Rotating Objective Turret with Power Indicator: has three lenses attached with varying powers of magnification.
 - The shortest lens (5x) has the weakest magnifying power.
 - The longest lens (60x) has the highest level of magnifying power.
 - If the 60x lens is used, and the zoom eyepiece is turned to maximum power, the object is under the microscope will be magnified 1200 times. It will appear 1200 times larger than when viewed with the naked eye.
- E. Stage: is located just below the rotating objective turret. The stage is a flat surface with a hole in the center to allow light to reflect from the mirror, or the battery-operated viewing light, to the microscope lenses.
- F. Stage Clips: hold the slide firmly in place.
- G. Reflecting Mirror / Viewing Light: are located just below the stage. The angle of the mirror can be adjusted to control the amount of light entering the microscope. The viewing light is located beneath the mirror. If additional light is needed, flip the mirror over and the light will turn on. When facing down, the light will be off.

Important: To switch from mirror to light, and vice versa, flip the light up and down. Do not turn the mirror/light from side to side in a clockwise or counterclockwise manner. If you inadvertently do so, the light will not go in the upright position. To correct, make a ½ rotation clockwise.

H. Horseshoe-Shaped Base: features a rubberized, non-skid bottom. Remove the base to replace the two AA batteries.

I. Projector Hood: projects and image of the specimen viewed on the stage into a wall or screen. Ideal for group study. (When not in use, unscrew the projector hood by turning counterclockwise. Replace it with the zoom eyepiece.

Battery installation: The light on the microscope requires fresh AA alkaline batteries. The battery compartment is located in the base of the unit. Carefully remove the cover on the base to insert or replace two AA batteries inside as shown in the diagram. Do not mix old and new batteries. Be sure to put the batteries in correctly.

USING THE PROJECTOR HOOD

The projector hood lets small or large groups of students view a specimen at the same time. The projector hood can be used to project an image in two ways: viewing screen on the hood itself (for small groups), and via wall or screen projection (for class wide viewing). Here is how to use the projector hood:

To attach the projector hood:

1. Remove the eyepiece by turning it counterclockwise.
2. Place the hood on top of the body tube.
3. Turn the projector hood clockwise until it is securely in place.

Using the projector hood with a small group

1. When using the projector hood, the light bulb, rather than the mirror, will be used. Flip the mirror down so the light bulb is upright projecting light through the stage.
2. Place a slide onto the microscope stage.
3. You will see the specimen projected on the viewing screen. Adjust the focus knobs until the image is sharp and clear.

Using the projector hood for projecting on a wall or screen

1. Attach a piece of paper to the wall at the same height as the projector hood.

2. Place the microscope 15 to 20 inches from the wall. Aim the viewing port toward the paper (or screen).
3. Remove the viewing screen from the projector hood.
4. Make the room as dark as possible. An image of the specimen on the stage will be projected on the wall.

INTRODUCTORY MICROSCOPE ACTIVITY

This activity will help students understand the process of the increasing magnification.

1. Remove the dust cap from the eyepiece.
2. Rotate the focus knob to raise the body tub to the highest point. Turn the turret to the lowest power of magnification (50x/100x).
3. Cut a single alphabet letter out of an old newspaper or magazine.
4. Place its letter on a clear slide.
5. Fill the dropper with water and place a drop on the slide.
6. Place a cover slip over the water and letter.
7. Using the stage clip, place the slide on the stage and move the "letter" towards the center holed.
8. Rotate the focus knob and lower the body tube to the lowest point. NOTE: TO AVOID BREAKAGE, DO NOT TOUCH THE SLIDE WITH THE VIEWING LENS.
9. Now rotate the focus knob slowly to raise the body tube until a clear image of the letter is seen through the eyepiece.
10. Adjust the angle of the mirror to let an adequate amount of light enter the microscope. (Use the light if necessary.) The image you will see is upside down and reversed. When you move the slide to the left, the image will move to the right and vice versa.
11. By rotating the objective turret, you can change the magnification of the image. Notice that when the image is enlarged, less of the subject and its illumination will decrease accordingly. Doubling the magnifying power reduces the brightness by 75%. Adjust the mirror to allow more light to enter the microscope.
12. Repeat its process of magnification several times until you are familiar with the process of increasing magnification.

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